

INNOVATION POWERED BY NATURE

SOLIS BIODYNE OÜ

E: info@solisbiodyne.com


T: +372 7409 960

F: +372 7402 079



Unique enzyme with
Stability TAG technology for increased temperature tolerance

solisbiodyne.com



[1] Kahre, O. et al., Compositions for increasing polypeptide stability and activity, and related methods, EP2501716B1 (2015) and US9321999B2 (2016).

[2] Jumper, J. et al. Highly accurate protein structure prediction with AlphaFold. Nature (2021).

[3] Varadi, M. et al. AlphaFold Protein Structure Database: massively expanding the structural coverage of protein-sequence space with high-accuracy models. Nucleic Acids Research (2021).

Features

- Stable at 25°C for at least 2 months and at 37°C 2 weeks
- Heat-labile (compatible with Sanger sequencing). No reactivation is detected after heat inactivation
- Fast 30 sec reaction time
- Tolerant to common inhibitors
- Reaction set-up and shipment without ice
- Glycerol-free formulation is available

Applications

- Widely used to eliminate carryover contamination in PCR and LAMP
- Enhancer of cloning efficiency of PCR products
- Site-directed mutagenesis
- As a probe for protein-DNA interaction studies
- Glycosylase-mediated single nucleotide polymorphism detection (GMPD)
- Study of DNA repair and mutation detection
- SNP genotyping




Figure 1. 3D model of Salini UNG™ Uracil-N-Glycosylase protein structure with the Stability TAG. Protein structure was predicted with AlphaFold 2.




Figure 2. Heat inactivation of Salini UNG™ Uracil-N-Glycosylase. UNG activity is measured at 37°C for 40 minutes using Bio-Rad CFX96 platform by the release of fluorescence from a uracil containing probe labeled with a FAM fluorophore and a quencher. In a resting state the probe forms a duplex and FAM fluorescence is quenched. An active UNG cleaves uracil, the duplex dissociates and FAM fluorescence is emitted. To test the inactivation of

the enzyme, Salini UNG™ was heat-treated at 70°C for 5 minutes. No reactivation was detected after storing the heat-treated samples at 4°C (Figure 2) or 25°C (data not shown) for 48h.

Prevention of carryover contamination



- Substitute dTTP with dUTP in your PCR or LAMP reaction. Preferably choose polymerases that can tolerate 100% substitution such as **Solis BioDyne SolisFAST® DNA Polymerase** (included in SolisFAST® (q)PCR mixes).
- Add Salini UNG™ to master mix.
- If you set up the reaction on ice add 30 sec at 25°C treatment step before amplification. No additional step is required if you set up reaction at room temperature.
- Salini UNG™ removes uracil (U) incorporated into single- or double-stranded DNA. Degraded dU-DNA cannot be used as a template for further amplification.
- Enjoy peace of mind not worrying about possible carryover contamination from repeated amplification of the same target sequence.






Figure 3. qPCR reactions were performed using SolisFAST® Probe qPCR Mix (no ROX) without UNG on the Bio-Rad CFX96 platform, and four 10-fold serial dilutions of human gDNA (from 2 ng/µl to 0.002 ng/µl) in green and NTC in red. For contamination simulation, approximately 4500 cp/ rxn of previously amplified and purified dU-DNA was used. Salini UNG™ Uracil-N-Glycosylase was added to the reaction (b) at a final concentration of 0.027 U/µl. Contaminated samples can result in false positives that compromise the data (a). Salini UNG™ Uracil-N-Glycosylase degrades 100% of contaminant dUTP-DNA (b) and does not compromise performance compared to the reaction with no contamination simulation and no UNG present in the mix (c).

Ordering information

Bulk solutions available

Application	Product	Size¹	CAT. NO.
Variable	Salini UNG™ Uracil-N-Glycosylase, 1 U/µl	25 µl 100 µl 5 x 100 µl	31-01-00000-S (free sample) 31-01-00100 31-01-00100-5
PCR	SolisFAST® Master Mix with UNG	250 rxn 5 x 250 rxn 5000 rxn	24-21-00001 24-21-00001-5 24-21-00020
	SolisFAST® Master Mix with UNG, Ready To Load	250 rxn 5 x 250 rxn 5000 rxn	24-22-00001 24-22-00001-5 24-22-00020
qPCR	SolisFAST® Probe qPCR Mix with UNG (ROX)	250 rxn 5 x 250 rxn 5000 rxn	28-21-00001 28-21-00001-5 28-21-00020
	SolisFAST® Probe qPCR Mix with UNG (no ROX)	250 rxn 5 x 250 rxn 5000 rxn	28-22-00001 28-22-00001-5 28-22-00020
	SolisFAST® Probe qPCR Mix with UNG (Purple)	250 rxn 5 x 250 rxn 5000 rxn	28-23-00001 28-23-00001-5 28-23-00020
RT-PCR	One-step RT-PCR Kit with UNG		Available upon request
RT-qPCR	SOLIScript® Fast 1-step RT-qPCR Mix with UNG		Launch in January 2023

 1 Volume in μl or count of 20 μl reactions. All PCR and qPCR Master Mixes are 5x-concentrated solutions supplied in 1 ml or 20 ml vials or bottles (a' 250 rxn and 5000 rxn, respectively).

FL-31-01-V1



For further details and ordering please contact info@solisbiodyne.com or call +372 740 9960

