

## SoliSD™ Lyo-compatible RT-LAMP Kit



Store at **-20°C**  
upon receipt

Cat. No.	Size (25 µl reactions)
32-23-0000S	250 reactions
32-23-01000	1000 reactions

### Description:

Isothermal DNA polymerase is an *in silico* designed DNA-directed DNA polymerase, based on the Large Fragment of the DNA polymerase from the genus *Bacillus smithii*. The enzyme has a strong strand displacement DNA polymerase activity and lacks 5'→3' and 3'→5' exonuclease activities. Thus, the enzyme is useful for applications where synthesis through double-stranded DNA regions is required, such as loop-mediated isothermal amplification (LAMP). Isothermal DNA Polymerase functional similarity to Bst and Bsm DNA polymerases allows its use in respective applications. Isothermal DNA polymerase is expressed recombinantly in *E. coli* and has enhanced stability due to the incorporated Stability TAG (covered by the patent EP2501716). SoliSD™ Glycerol-Free Bsm DNA Polymerase (40 U/µl) is fully compatible with lyophilization.

### Reagents supplied:

Component	Catalogue Number	
	32-23-0000S	32-23-01000
SoliSD™ Bsm DNA polymerase (40 U/µl)	50 µl	200 µl
10x RT-LAMP Reaction Buffer 200 mM Tris-HCl; 150 mM (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> ; 1% Tween-20; 250 mM KCl; additives	1 ml	3 ml
300x RT mix containing FIREScript®, SOLIScript®, RiboGrip®	250 rxn	1000 rxn
100 mM MgSO <sub>4</sub>	500 µl	2 ml
10x GC-rich Enhancer	1 ml	3 ml

### Shipping and Storage conditions:

Routine storage: at -20°C.

The recommended reaction mixture protocol and cycling conditions are based on Solis BioDyne test results and should be considered as a starting point for reaction optimization tests. The optimal conditions may vary depending on the cycler model and ramp rate.

Manufactured by Solis BioDyne in compliance with the ISO 9001 and ISO 13485 certified Quality Management System.

**FOR RESEARCH USE ONLY.  
NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

### Recommended RT-LAMP reaction mix with SoliSD™ Bsm DNA Polymerase:

Component	Volume	Final concentration
Water, nuclease-free	up to 25 µl	
10x RT-LAMP Reaction Buffer	2.5 µl	1x
100 mM MgSO <sub>4</sub> <sup>a</sup>	1.6 µl	6.4 mM
10x GC-rich Enhancer <sup>b</sup>	Variable	Variable
SoliSD™ Glycerol-Free Bsm DNA Polymerase (40 U/µl) <sup>c</sup>	0.2 µl	0.32 U/µl
300x RT Mix <sup>d</sup>	0.083 µl	1x
dNTPs	Variable	1-2 mM
RT-LAMP primers <sup>e</sup>	Variable	Variable
RNA template	Variable	Variable
<b>Total</b>	<b>25 µl</b>	

<sup>a</sup> Optimal concentration of MgSO<sub>4</sub> in LAMP is 6-7 mM, consider higher concentrations for multiplex reactions.

<sup>b</sup> The addition of GC-rich Enhancer is optional, consider adding 1x, 2x, or 5x final conc. to the reaction mix in case of a GC-rich template.

<sup>c</sup> The optimal concentration of polymerase should be titrated depending on the application.

<sup>d</sup> Lyo-compatible reverse transcriptase mix with RNase Inhibitor.

<sup>e</sup> Optimal results may require titration of primer concentration, with a starting point of 200 nM for F/B, 400 nM for loop primers, and 1,6 µM for inner primers.

For real-time fluorescence-based product detection, add intercalating fluorescent dye into the reaction mix according to the manufacturer's recommendation.

### Perform reaction at 59°C for 20-60 min.

SoliSD™ Bsm DNA polymerase working temperature is 51-62°C. The optimal temperature of the assay may vary depending on the primers and final buffer conditions.

We recommend screening several sets of primers to achieve the best sensitivity and eliminate the possibility of amplification in an NTC (no-template control) reaction.

### Contamination control:

Due to the nature of the enzyme and primers used in LAMP reactions, opening reaction vessels following amplification is not recommended.

**Safety warnings and precautions:**

This product and its components should be handled only by persons trained in laboratory techniques. It is advisable to wear suitable protective clothing, such as laboratory overalls, gloves, and safety glasses. Care should be taken to avoid contact with skin or eyes. In case of contact with skin or eyes, wash immediately with water. Refer to the Safety Data Sheet for more information.

**Technical support:**

Contact your sales representative for any questions or send an email to [support@solisbiodyne.com](mailto:support@solisbiodyne.com)

Online chat is available at [www.solisbiodyne.com](http://www.solisbiodyne.com)

DS-32-23 v2

Revision date: 17.09.2024

**Permitted Use:** This product is supplied for Research Use Only (the **Permitted Use**). It is suitable for use as a component of molecular diagnostic assays, where applicable country laws allow. This product alone does not provide any diagnostic result. If the customer wishes to use the product for any purpose other than the Permitted Use, including (without limitation) resale or alteration, the customer should obtain the appropriate license from Solis BioDyne. Some applications of this product may require a license/licenses from one or more third parties which are not provided by the purchase of this product. Users should obtain the license if required. Covered by the patent EP2501716, made by the methods of US Patent No 9.321.999.

**Trademark information:** SoliSD™ is a trademark of Solis BioDyne OÜ. FIREScript®, SOLIScript® and RiboGrip® are registered trademarks of Solis BioDyne OÜ.

**Warranty and Disclaimer:** This product shall comply with its relevant specification and be fit for its stated purpose, but Solis BioDyne gives no other warranty and makes no representation as to description or quality. Any such warranty or representation is excluded, to the fullest extent permitted by law. In particular, but without limiting the foregoing, Solis BioDyne shall not be liable for the failure of the product to comply with its relevant specification where such failure arises as a result of: (i) customer negligence or because the customer failed to follow any of the applicable technical data or safety sheets, standard user materials, use guidelines or any other information provided by Solis BioDyne as to the storage, transportation, handling, use or maintenance of the products or other good practice regarding the same, or (ii) the customer altering the products in any way without the prior written consent from Solis BioDyne, or (iii) the products differing from the relevant specification as a result of changes made to ensure their compliance with applicable statutory or regulatory requirements.

Nothing shall limit or exclude Solis BioDyne's liability for death or personal injury caused by its negligence, fraud or fraudulent misrepresentation or any matter in respect of which it would be unlawful for Solis BioDyne to exclude or restrict liability. Without limiting the foregoing, Solis BioDyne shall under no circumstances whatever be liable to the customer, whether in contract, tort (including negligence), breach of statutory duty, or otherwise, for any loss of profit, or any indirect or consequential loss arising under or in connection with the products and Solis BioDyne's total liability to the customer in respect of all other losses arising under or in connection with the product, whether in contract, tort (including negligence), breach of statutory duty, or otherwise, shall in no circumstances exceed the price of the products supplied in respect of which the liability has arisen.

**Solis BioDyne**

Teaduspargi 9, 50411 Tartu, Estonia, tel: +372 740 9960, fax: +372 740 2079, e-mail: [info@solisbiodyne.com](mailto:info@solisbiodyne.com), [www.solisbiodyne.com](http://www.solisbiodyne.com)