

SoliSD™ Lyo-compatible *Bsm* DNA polymerase Kit

Catalogue Number	Size
32-22-0000S	100 reactions
32-22-00250	250 reactions
32-22-01000	1000 reactions



**Store at –20°C
upon receipt**

Shipping:

At room temperature

Batch Number and Expiry Date:

See vial

Storage and Stability*:

- Routine storage at –20°C (–28°C to –18°C) until Expiry Date
- Stable at 4°C (2°C to 8°C) for up to 6 months
- Stable at room temperature (15–25°C) for at least 1 month
- Freeze-thaw stability: 10 cycles

Reaction setup:

At room temperature

Manufactured by Solis BioDyne, in compliance with the ISO 9001 and ISO 13485 certified Quality Management System.

Product description:

- SoliSD™ *Bsm* DNA polymerase is based on the Large Fragment of the DNA polymerase from the genus *Bacillus smithii* with functional similarity to other *Bst* and *Bsm* DNA polymerases. It has enhanced stability due to the incorporated Stability TAG (covered by the patent EP2501716)**.
- The enzyme has a strong strand displacement DNA polymerase activity and lacks 5'→3' and 3'→5' exonuclease activities. It is designed for applications requiring synthesis through double-stranded DNA regions, such as loop-mediated isothermal amplification (LAMP).
- Unique SoliSD™ Supplement system:
 - Resolves common NTC signal issue
 - Minimizes variability between replicates
 - Allows reaction set-up at room temperature
- SoliSD™ *Bsm* DNA polymerase Kit is provided in a flexible 5-vial format. The development of an isothermal amplification assay can be challenging, and optimization depending on the target, primer set, and detection method is often required.

Features

1. Enzyme is active at a wide temperature range between 51-62°C with optimum at 60°C
2. No 5'→3' and 3'→5' exonuclease activity
3. Strong strand displacement DNA polymerase activity
4. Short 4–20-minute reaction time to result
5. No NTC signal

6. Increased thermostability for at least 1 month at 37°C

7. Market-level inhibitor tolerance

Reagents supplied:

Component	Catalogue Number		
	32-22-0000S	32-22-00250	32-22-01000
SoliSD™ <i>Bsm</i> DNA polymerase (40 U/μl)	20 μl	50 μl	200 μl
10x Isothermal Reaction Buffer (pH 8.8); 200 mM Tris-HCl, 150 mM (NH ₄) ₂ SO ₄ , 1% Tween-20, 250 mM KCl	1 ml	1 ml	3 ml
100 mM MgSO ₄	500 μl	500 μl	2 ml
25x SoliSD™ Supplement	250 μl	250 μl	1 ml
10x GC-rich Enhancer	500 μl	1 ml	3 ml

Step-by-step guidelines for wet use:

1. Thaw the reagents at room temperature. Mix by vortexing or pipetting up and down, then centrifuge briefly.
2. Prepare a reaction master mix. Add components except for template DNA in the suggested order. Mix by vortexing or pipetting up and down, then centrifuge briefly.

Reaction component	Volume	Final conc.
Nuclease-free water ^a	Up to 25 μl	
10x Isothermal Reaction Buffer	2.5 μl	1x
100 mM MgSO ₄	1.6 μl	6.4 mM ^c
10x GC-rich Enhancer ^b	optional	optional

dNTPs ^a	variable	1.4 mM ^d
25x SoliSD™ Supplement	1 µl	1x
Fluorescent Dye ^a	variable	0.3-0.5 µM
FIP/BIP primers ^a	variable	1.6 µM
F3/B3 primers ^a	variable	0.2 µM
LoopF/B primers ^a	variable	0.4 µM
SoliSD™ <i>Bsm</i> DNA polymerase (40 U/µl)	0.2 µl	0.32 U/µl
Template DNA ^a	variable	10 ¹ -10 ¹⁰ cp/rxn
Total	25 µl	

^a Not included in the kit; ^b Addition of GC-rich Enhancer is optional, consider adding 1x, 2x, or 5x final concentration to the reaction mix in case of a GC-rich template.

^c Optimal concentration of MgSO₄ in LAMP is between 6-7 mM. ^d Optimal concentration of dNTPs is between 1-2 mM.

3. Dispense the reaction mixture and add DNA template.
4. Perform the reaction at 60°C for 10-45 min. SoliSD™ *Bsm* DNA polymerase is active at 51-62°C. The optimal temperature of the assay may vary depending on the primers and final buffer conditions.
5. If need be, inactivate the enzyme by heating it at 80°C for 10 min.

Step-by-step guidelines for lyophilization:

1. Thaw the reagents at room temperature. Mix by vortexing or pipetting up and down, then centrifuge briefly.
2. Prepare a lyophilization master mix. **Add components except for 10x Isothermal Reaction Buffer, 100 mM MgSO₄, 10x GC-rich Enhancer, and template DNA in the suggested order.**

Reaction component	Volume	Final conc.
Nuclease-free water ^a	Up to 25 µl	
dNTPs ^a	variable	1.4 mM
25x SoliSD™ Supplement	1 µl	1x
Fluorescent Dye ^a	variable	0.3-0.5 µM
FIP/BIP primers ^a	variable	1.6 µM
F3/B3 primers ^a	variable	0.2 µM
LoopF/B primers ^a	variable	0.4 µM
Excipient mix ^a	variable	1x
SoliSD™ <i>Bsm</i> DNA polymerase (40 U/µl)	0.2 µl	0.32 U/µl
Total	25 µl	

^a Not included in the kit.

3. Vortex and centrifuge the reaction mix briefly. Dispense appropriate volumes of the reaction mix to lyophilization vials or other suitable containers (e.g., PCR tubes, plates).
4. Proceed with lyophilization immediately. You can find the recommended lyophilization protocol on the product webpage under Documents section or proceed based on the equipment you are using.

Step-by-step guidelines for using lyophilized product:

1. Prepare 2x reconstitution mix using 10x Isothermal Reaction Buffer, 100 mM MgSO_4 , and 10x GC-rich Enhancer (optional) (it is recommended to optimize MgSO_4 and 10x GC-rich Enhancer concentrations empirically). The mix may also be prepared at concentrations other than 2x, but care must be taken to ensure that the final concentration in a 25 μl reaction are 1x for 10x Isothermal Reaction Buffer and 6-7 mM for MgSO_4 after addition of both reconstitution mix and template DNA.

Reaction component	Volume	Final conc.
10x Isothermal Reaction Buffer	2.5 μl	2x
100 mM MgSO_4^{a}	1.6 μl	12.8 mM
10x GC-rich Enhancer ^b	Optional (0-8.4 μl)	Optional
Nuclease-free water	Up to 12.5 μl	

^a Optimal concentration of MgSO_4 in LAMP is between 6-7 mM. ^b Addition of GC-rich Enhancer is optional, consider in case of a GC-rich template; final concentration of 1-3x in 25 μl is usually used.

2. Reconstitute the lyophilized cake or bead by adding 12.5 μl of reconstitution mix.
3. Vortex until cake or bead is solubilized, then spin down.
4. Add 12.5 μl of DNA template/nuclease-free water up to make a total reaction volume of 25 μl .
5. Perform the reaction at 60°C for 10-45 min. SoliSD™ *Bsm* DNA polymerase is active at 51-62°C. The optimal temperature of the assay may vary depending on the primers and final buffer conditions.
6. If need be, inactivate the enzyme by heating it at 80°C for 10 min.

Primer design and validation

LAMP assay requires a set of 4-6 primers. We recommend screening several sets of primers to achieve the best sensitivity and eliminate the possibility of the signal in an NTC (no-template control) reaction. To facilitate primer design use dedicated primer design software.

Contamination control

To prevent contamination, we recommend the separation of working areas into pre-amplification (reaction mixture preparation), amplification (dispensing of reaction mixture and addition of template), and post-amplification. Note! Due to the nature of enzyme and primers used in LAMP reactions, opening reaction vessels following amplification is not recommended.

Controls

We recommend including several controls:

- Negative controls: no template and non-target template controls.
- Positive control: a previously validated template and primer set.

Source:

Purified from an *E. coli* strain that carries an overproducing plasmid containing a SoliSD™ *Bsm* DNA polymerase gene.

Unit definition:

The unit is determined in DNA strand displacement assay, where 1 unit of the enzyme catalyzes displacement of 40 nmol of oligonucleotide probe from the complementary DNA strand in 30 min at 60 °C.

Safety precautions:

Please refer to Safety Data Sheet for more information.

Technical support:

Contact your sales representative for any questions or send an email to support@solisbiodyne.com

For research use only. Not for use in diagnostic procedures.

DS-32-22 v4. Revised 07.11.2025

***Product stability** is assessed using routine QC assays and QC criteria set forth in the product specification and are intended to provide guidelines for shipping and storage conditions only. Customer or its designee shall be responsible for conducting all necessary stability testing applicable to their assay and/or QC criteria, and to comply with any applicable regulatory requirements or guidelines. Such stability testing shall include testing to validate the lead times for shipment, the shelf life of, and the product specifications applicable to shipment, storage and handling of the assay assembled and packed by the customer.

Permitted Use: This product is supplied for research use only (the Permitted Use). Some applications of this product may require a license/licenses from Solis BioDyne OÜ or one or more third parties which are not provided by the purchase of this product. This product shall comply with its relevant specification and be fit for its stated purpose, Solis BioDyne OÜ gives no other warranty and makes no representation as to description or quality. For more information and full disclaimers, please contact our customer service.

****Covered by patent EP2501716, made following the methods of US Patent No 9,321,999.**

Trademark information: SoliSD™ is a trademark of Solis BioDyne OÜ.

Manufacturer: Solis BioDyne OÜ | Teaduspargi 9, 50411 | Tartu, Estonia (EU)