

SOLIS BIODYNE OÜ

E: info@solisbiodyne.com T: +372 7409 960 F: +372 7402 079

SolisFAST® Master Mix

Fast singleplex and multiplex endpoint <u>PCR</u> assays!

solisbiodyne.com



## **Features**

- Fast less time from sample to results
- Accurate reliable results in single- and multiplex assays
- Stable 30-day stability at room temperature
- Trustworthy avoid carry-over contamination with UNG mixes and get reliable results in Sanger sequencing applications

The mix contains all components necessary, except primers, nuclease-free water and DNA template.

SolisFAST® Master Mix is a 5x-concentrated, ready-to-use solution for fast singleplex and multiplex endpoint PCR assays. The mix allows amplification of up to 5 kb fragments from low complexity DNA templates (e.g. lambda, plasmid DNA), and up to 3 kb from complex DNA templates (e.g. human genomic DNA). Combining our novel in-silico designed SolisFAST® DNA Polymerase and optimized buffer, this mix offers robust PCR, accurate target detection, ice-free shipping and reaction set-up.

#### Results in 20 minutes!




Figure 1. SolisFAST® Master Mixes are suitable for 'slow' and 'fast' PCR cyclers. HIRA gene fragment (515 bp) from human gDNA was amplified using SolisFAST® Master Mix (lane 1-2) and Master Mix with UNG (lane 3-4). Fast cycling settings (initial activation 1 min at 98 °C; denaturation 5 sec at 98 °C, annealing/extension 20 sec at 60 °C, 25 cycles) were used on Biometra T1 Thermocycler (lane 1 and 3) and Eppendorf® Mastercycler® X50s (lane 2 and 4).


PCR cycler	Cycler's ramp rate	PCR run time
Biometra T1 = 'slow' machine	4 °C/sec	26 min
Eppendorf X50s = 'fast' machine	10 °C/sec	20 min

Worldwide ice-free shipping!

# Excellent singleplex and multiplex PCR results!

Figure 2. Equal performance of SolisFAST® Master Mix in singleplex and multiplex assay.

Eighteen fragments of human gDNA (ranging from 122 to 1340 bp) were amplified in singleplex and multiplex PCR. Amplifications were performed as a single run on Eppendorf® Mastercycler® X50s using 3-step cycling program optimized for multiplex assays (initial activation 2 min at 98 °C; denaturation 10 sec at 98 °C, annealing 10 sec at 60 °C, extension 30 sec at 72 °C (30 cycles).



### dUTP and UNG don't affect Sanger sequencing results!



Figure 3. Sanger sequencing of products obtained with SolisFAST® Master Mix and Master Mix with UNG. The presence of UNG enzyme in PCR sample doesn't affect the sequencing results after 24-hour incubation of PCR products at room temperature.

Avoid carry-over contamination with UNG Mixes!

### **Ordering information**

Product	CAT. NO.	RXN / 20 μl
SolisFAST® Master Mix	24-01-0000S (free sample) 24-01-00001 24-01-00001-5 24-01-00020	50 rxn / 20 μl 250 rxn / 20 μl 5x (250 rxn / 20 μl) 5000 rxn / 20 μl
SolisFAST® Master Mix, Ready To Load	24-02-0000S (free sample) 24-02-00001 24-02-00001-5 24-02-00020	50 rxn / 20 μl 250 rxn / 20 μl 5x (250 rxn / 20 μl) 5000 rxn / 20 μl
SolisFAST® Master Mix with UNG	24-21-0000S (free sample) 24-21-00001 24-21-00001-5 24-21-00020	50 rxn / 20 μl 250 rxn / 20 μl 5x (250 rxn / 20 μl) 5000 rxn / 20 μl
SolisFAST® Master Mix with UNG, Ready To Load	24-22-0000S (free sample) 24-22-00001 24-22-00001-5 24-22-00020	50 rxn / 20 μl 250 rxn / 20 μl 5x (250 rxn / 20 μl) 5000 rxn / 20 μl



